

# Incline Platform Lifts & Stair Glides



## IL500 Sierra Incline Platform Lift:

The innovative Sierra IL500 is ideal for transporting a user and their manual wheelchair up a straight stairway. Can use most wheelchairs up to 36" in length and offers an integrated seat option for use as a stair lift.

### Features & Benefits:

- 500-lb capacity
- Accepts most wheelchairs up to 36" in length
- Space saving: 13" folded when not in use
- Free Standing Post Kit for NO wall modifications
- Safety Sensors stop lift in case of an obstacle
- Easy Installation with no load-bearing wall needed.
- Residential Use Only

### Optional Accessories:

- Integrated seat option for use as a stair lift
- Multiple platform configurations
  - 27.5" X 36"
  - 27.5" X 36" 90-degree exit platform

### Safety Features:

- Key lock
- Emergency stop button
- Additional foot of track (Track, Handrail, & Mounting Brackets)

## BMS IS THE MOBILITY EXPERT!

**We custom build to your specifications.**

**Call BMS FIRST for a FREE No-Obligation Personal Mobility Evaluation!**  
**1-800-443-3390**

**Boardman Medical Supply**  
**"Visit one of our Locations in OH & PA"**  
**Toll Free: 800-443-3390**  
**Local: 330-545-6700**  
**Fax: 800-952-5352**

GIRARD, OH (Corporate) \*  
 300 N. State Street  
 Girard, OH. 44420  
 Phone: (330) 545-6700

EAST LIVERPOOL, OH  
 500 Market St.  
 East Liverpool, OH 43920  
 Phone: (330) 385-6412  
 (By Appointment Only)

BOARDMAN, OH \*  
 404 Boardman-Canfield Rd.  
 Boardman, OH 44512  
 Phone: (330) 758-0595

BEACHWOOD, OH  
 24000 Merchantile Rd, Ste 6  
 Beachwood, OH 44122  
 Phone: (330) 998-6012  
 (By Appointment Only)

CANTON, OH  
 1204 30th St. N.W.  
 Canton, OH 44709  
 Phone: (330) 492-4300

WARREN, OH  
 1309 E. Market St.  
 Warren, OH 44483  
 Phone: (330) 395-7252  
 (By Appointment Only)

Dedicated to Quality Service & Equipment



ACHC ACCREDITED



YOUR CHOICE . . . OUR PROMISE  
 800-443-3390

[www.boardmanmedicalsupply.com](http://www.boardmanmedicalsupply.com)

Fax: 800-952-5352

## HOME STAIR GLIDES VEHICLE MOBILITY LIFTS PORCH LIFTS & RAMPS



## BMS IS THE MOBILITY EXPERT!

It's easy and affordable to take your mobility product with you wherever you go.

**Call BMS FIRST for a FREE Personal Mobility Evaluation!**

BOARDMAN • MEDICAL • SUPPLY



COMPANY, INC.

YOUR CHOICE . . . OUR PROMISE

800-443-3390

[www.boardmanmedicalsupply.com](http://www.boardmanmedicalsupply.com)

Fax: 800-952-5352

# Vehicle Lifts

## Outside Lifts



- Carries virtually all power chairs
- 4 – Self-tensioning, Q'Straint retractors
- No power chair modifications
- Manual back-up
- License plate mount
- Swing-Away option available
- Made in the USA

## Inside Lifts



- Carries virtually all standard power chairs and scooters-- no mobility device modifications required
- Most compact unit on the market-- requires just 43" of depth
- Allows continued use of second row seating; no drilling required in most vehicles
- Manual crank backup
- Weatherproof, remote hand control (wireless option available)
- Ideal for most minivans, SUVs, and crossovers
- Made in the USA



# Platforms / Ramps

## Vertical Platforms



- Automatic, folding access ramp
- Available in AC and DC models
- Solid platform with non-skid surface for added safety
- Solid 36" high guard panels
- Weather-resistant design for indoor and outdoor use
- Constant pressure, paddle switch controls with emergency stop button
- Sensitive safety pan on bottom of platform stops the lift if it meets an obstruction
- Made in USA

## Modular Ramps



- Multiple modular configurations
- Non-skid surface safety feature
- Fits any unique entryway situation
- Easy install...no permits necessary
- Maintenance free aluminum
- Gates available
- Will not rot or rust

## Portable Ramps



- Portable
- Durable
- Inexpensive
- Use with scooters and wheelchairs

# Lifts / Glides

## Stair Glides



- Proven Rack and Pinion drive system
- Quiet Drive System with Soft Start and Stop
- DC Drive even runs during power interruptions
- Aesthetic track, finished floor mounts, and integrated charging bases
- Upholstered seat and armrest w/ folding footrest
- Includes two wireless remote controls
- Two-piece chassis cover design
- Easy access to manual back-up
- 300 lbs. capacity

## Pool Lifts



- 350 lb. capacity
- Meets and exceeds ADA specifications
- Ambidextrous seating
- Secure concrete anchor mounting with sliding tube for removal and storage
- Fully submersible, comfortable seating
- Durable stainless steel construction
- Folding arms & foot plate
- Sealed hand control
- Wall-mounted charger with hanger
- Removable, easy-access battery
- Optional pool lift cover available